

Product Datasheet

NGFR peptide (orb217675)

Description	NGFR peptide
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	NGFR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214319.
Purity	> 85%
Source	Synthetic
Uniprot ID	P07174, Q9Z0W1, P08138
Entrez	24596, 4804
Expiration Date	12 months from date of receipt.