

Product Datasheet

NSF peptide (orb217693)

Description	The peptide is used to block anti-NSF antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	NSF
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214337.
Purity	> 85%
Source	Synthetic
Uniprot ID	P46459, Q9QUL6, P46460
Entrez	60355, 18195, 4905
Expiration Date	12 months from date of receipt.