

Product Datasheet

TrkB peptide (orb217695)

Description	TrkB peptide
Reactivity	Human, Mouse, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	NTRK2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214339.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15209, Q63604, Q16620
Entrez	4915, 18212, 25054
Expiration Date	12 months from date of receipt.