

Product Datasheet

OAT peptide (orb217699)

Description	The peptide is used to block anti-OAT antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	OAT
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214343.
Purity	> 85%
Source	Synthetic
Uniprot ID	P04181, P29758, P04182
Entrez	4942, 18242, 64313
Expiration Date	12 months from date of receipt.