

Product Datasheet

PC peptide (orb217720)

Description	The peptide is used to block anti-Pyruvate Carboxylase antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	PC
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214364.
Purity	> 85%
Source	Synthetic
Uniprot ID	P11498, P52873, Q05920
Entrez	5091, 25104
Expiration Date	12 months from date of receipt.