

Product Datasheet

PCNA peptide (orb217723)

Description	The peptide is used to block anti-PCNA antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	PCNA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214367.
Purity	> 85%
Source	Synthetic
Uniprot ID	P04961, P12004, P17918
Entrez	5111, 18538, 25737
Expiration Date	12 months from date of receipt.