

Product Datasheet

PDE4B peptide (orb217726)

Description	The peptide is used to block anti-PDE4B/C/D antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	PDE4B; PDE4C; PDE4D
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214370.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q08493 , Q08499 , Q07343
Entrez	5144 , 5143 , 5142
Expiration Date	12 months from date of receipt.