

Product Datasheet

PDGFRB peptide (orb217730)

Description	The peptide is used to block anti-PDGFR beta antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	PDGFRB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214374.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05622, Q05030, P09619
Entrez	24629, 5159, 18596
Expiration Date	12 months from date of receipt.