

Product Datasheet

CD31 peptide (orb217735)

Description	CD31 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PECAM1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214379.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q08481, Q3SWT0, P16284
Entrez	5175, 18613, 29583
Expiration Date	12 months from date of receipt.