

Product Datasheet

PEX7 peptide (orb217736)

Description	The peptide is used to block anti-PEX7 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PEX7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214380.
Purity	> 85%
Source	Synthetic
Uniprot ID	O00628, P97865
Entrez	18634, 5191
Expiration Date	12 months from date of receipt.