

Product Datasheet

PFKFB1 peptide (orb217738)

Description	The peptide is used to block anti-PFKFB1/4 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat, Zebrafish
Tested Applications	Blocking
Target	PFKFB1; PFKFB4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214382.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6DTY7, Q16877, P07953, P25114, P70266, P16118
Entrez	270198, 5210, 54283, 5207, 18639, 24638
Expiration Date	12 months from date of receipt.