

Product Datasheet

PFKL peptide (orb217740)

Description	The peptide is used to block anti-PFKL antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PFKL
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214384.
Purity	> 85%
Source	Synthetic
Uniprot ID	P12382, P17858, P30835
Entrez	5211, 18641, 25741
Expiration Date	12 months from date of receipt.