

Product Datasheet

PFN2 peptide (orb217742)

Description	The peptide is used to block anti-Profilin 2 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	PFN2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214386.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9EPC6 , Q9JJV2 , P35080
Entrez	18645 , 100909840 , 5217 , 81531
Expiration Date	12 months from date of receipt.