

Product Datasheet

SERPINB5 peptide (orb217749)

Description	The peptide is used to block anti-Serpin B5 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SERPINB5
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214393.
Purity	> 85%
Source	Synthetic
Uniprot ID	P36952, P70124, P70564
Entrez	20724, 5268, 116589
Expiration Date	12 months from date of receipt.