

Product Datasheet

PIK3C2G peptide (orb217751)

Description	The peptide is used to block anti-PI3K C2 gamma antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	PIK3C2G
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214395.
Purity	> 85%
Source	Synthetic
Uniprot ID	070173, 070167, 075747
Entrez	18705, 116720, 5288
Expiration Date	12 months from date of receipt.