

Product Datasheet

PLCG1 peptide (orb217761)

Description	The peptide is used to block anti-PLC gamma 1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	PLCG1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214405.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10686, Q62077, P19174
Entrez	5335, 18803, 25738
Expiration Date	12 months from date of receipt.