

Product Datasheet

PLCG2 peptide (orb217763)

Description	The peptide is used to block anti-PLC gamma 2
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PLCG2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214407.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8CIH5, P24135, P16885
Entrez	234779, 5336, 29337
Expiration Date	12 months from date of receipt.