

Product Datasheet

POLD1 peptide (orb217771)

Description	The peptide is used to block anti-DNA Polymerase delta 1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	POLD1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214415.
Purity	> 85%
Source	Synthetic
Uniprot ID	P52431, P28340, O54747
Entrez	5424, 59294, 18971
Expiration Date	12 months from date of receipt.