

## Product Datasheet

### POU2F1 peptide (orb217775)

|                            |                                                             |
|----------------------------|-------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to block anti-OCT1 antibody reactivity. |
| <b>Reactivity</b>          | Human, Porcine                                              |
| <b>Tested Applications</b> | Blocking                                                    |
| <b>Target</b>              | POU2F1                                                      |
| <b>Form/Appearance</b>     | Lyophilized powder                                          |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                |
| <b>Note</b>                | For research use only                                       |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb214419.    |
| <b>Purity</b>              | > 85%                                                       |
| <b>Source</b>              | Synthetic                                                   |
| <b>Uniprot ID</b>          | <b>P14859</b>                                               |
| <b>Entrez</b>              | <b>5451</b>                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                             |