

Product Datasheet

POU2F2 peptide (orb217776)

Description	The peptide is used to block anti-OCT2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	POU2F2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214420.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q00196, P09086
Entrez	5452, 18987
Expiration Date	12 months from date of receipt.