

Product Datasheet

PRKACA peptide (orb217791)

Description	The peptide is used to block anti-PKA C alpha/beta/gamma antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	PRKACA; PRKACB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214435.
Purity	> 85%
Source	Synthetic
Uniprot ID	P68181 , P68182 , P22612 , P22694 , P05132 , P27791 , P17612
Entrez	5567 , 18749 , 293508 , 5566 , 18747
Expiration Date	12 months from date of receipt.