

Product Datasheet

PRKAG1 peptide (orb217793)

Description	PRKAG1 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PRKAG1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214437.
Purity	> 85%
Source	Synthetic
Uniprot ID	O54950, P80385, P54619
Entrez	5571, 25520, 19082
Expiration Date	12 months from date of receipt.