

Product Datasheet

PRKAR1B peptide (orb217795)

Description	The peptide is used to block anti-PRKAR1B antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	PRKAR1B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214439.
Purity	> 85%
Source	Synthetic
Uniprot ID	P81377, P31321, P12849
Entrez	5575, 19085
Expiration Date	12 months from date of receipt.