

Product Datasheet

PRKAR2B peptide (orb217796)

Description	The peptide is used to block anti-PRKAR2B antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	PRKAR2B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214440.
Purity	> 85%
Source	Synthetic
Uniprot ID	P31323 , P12369 , P31324
Entrez	24679 , 19088 , 5577
Expiration Date	12 months from date of receipt.