

Product Datasheet

PRKCD peptide (orb217798)

Description	The peptide is used to block anti-PKC delta
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	PRKCD
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214442.
Purity	> 85%
Source	Synthetic
Uniprot ID	P28867, Q05655, P09215
Entrez	5580, 18753, 170538
Expiration Date	12 months from date of receipt.