

Product Datasheet

ERK1/2 peptide (orb217808)

Description	ERK1/2 peptide
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	MAPK1; MAPK3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214452.
Purity	> 85%
Source	Synthetic
Uniprot ID	P63086, P21708, P27361, P28482
Entrez	50689, 116590, 5595, 5594
Expiration Date	12 months from date of receipt.