

Product Datasheet

MAPK10 peptide (orb217811)

Description	The peptide is used to block anti-JNK3 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MAPK10
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214455.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q61831, P49187, P53779
Entrez	5602, 26414, 25272
Expiration Date	12 months from date of receipt.