

Product Datasheet

MASP1 peptide (orb217820)

Description	The peptide is used to block anti-MASP1 HC antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	MASP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214464.
Purity	> 85%
Source	Synthetic
Uniprot ID	P48740, P98064, Q8CHN8
Entrez	64023, 17174, 5648
Expiration Date	12 months from date of receipt.