

Product Datasheet

PTGS2 peptide (orb217832)

Description	The peptide is used to block anti-PGHS-2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	PTGS2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214476.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35355, P35354, Q05769
Entrez	19225, 5743, 29527
Expiration Date	12 months from date of receipt.