

Product Datasheet

PTK2 peptide (orb217835)

Description	The peptide is used to block anti-Focal Adhesion Kinase
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	PTK2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214479.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q05397 , O35346 , P34152
Entrez	5747 , 14083 , 25614
Expiration Date	12 months from date of receipt.