

Product Datasheet

RAD23B peptide (orb217845)

Description	The peptide is used to block anti-RAD23B antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	RAD23B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214489.
Purity	> 85%
Source	Synthetic
Uniprot ID	P54727, Q4KMA2, P54728
Entrez	298012, 5887, 19359
Expiration Date	12 months from date of receipt.