

Product Datasheet

RALB peptide (orb217849)

Description	The peptide is used to block anti-RALB antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	RALB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214493.
Purity	> 85%
Source	Synthetic
Uniprot ID	P11234, P36860, Q9JIW9
Entrez	5899, 116546, 64143
Expiration Date	12 months from date of receipt.