

Product Datasheet

NFkB p65 peptide (orb217859)

Description	NFkB p65 peptide
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	RELA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214503.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q04207 , Q04206
Entrez	19697 , 5970
Expiration Date	12 months from date of receipt.