

Product Datasheet

RFC2 peptide (orb217869)

Description	The peptide is used to block anti-RFC2 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	RFC2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214513.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9WUK4, P35250, Q641W4
Entrez	116468, 19718, 5982
Expiration Date	12 months from date of receipt.