

Product Datasheet

RGS10 peptide (orb217873)

Description	The peptide is used to block anti-RGS10 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	RGS10
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214517.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9CQE5 , P49806 , O43665
Entrez	67865 , 6001
Expiration Date	12 months from date of receipt.