

Product Datasheet

ROCK1 peptide (orb217876)

Description	The peptide is used to block anti-ROCK1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	ROCK1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214520.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q63644, P70335, Q13464
Entrez	6093, 81762, 19877
Expiration Date	12 months from date of receipt.