

Product Datasheet

RORC peptide (orb217877)

Description	RORC peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	RORC
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214521.
Purity	> 85%
Source	Synthetic
Uniprot ID	P51449, P51450
Entrez	19885, 6097
Expiration Date	12 months from date of receipt.