

Product Datasheet

RPL17 peptide (orb217882)

Description	The peptide is used to block anti-RPL17 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	RPL17
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214526.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9CPR4, P24049, P18621
Entrez	6139, 291434, 100526842
Expiration Date	12 months from date of receipt.