

## Product Datasheet

### RPS6KB1 peptide (orb217897)

|                            |                                                             |
|----------------------------|-------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to block anti-S6K1 antibody reactivity. |
| <b>Reactivity</b>          | Human, Mouse, Porcine, Primate, Rabbit, Rat                 |
| <b>Tested Applications</b> | Blocking                                                    |
| <b>Target</b>              | RPS6KB1                                                     |
| <b>Form/Appearance</b>     | Lyophilized powder                                          |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                |
| <b>Note</b>                | For research use only                                       |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb214541.    |
| <b>Purity</b>              | > 85%                                                       |
| <b>Source</b>              | Synthetic                                                   |
| <b>Uniprot ID</b>          | <b>P23443, Q8BSK8, P67999</b>                               |
| <b>Entrez</b>              | <b>6198, 72508, 83840</b>                                   |
| <b>Expiration Date</b>     | 12 months from date of receipt.                             |