

Product Datasheet

RPS6KB1 peptide (orb217898)

Description	The peptide is used to block anti-S6K1 antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat
Tested Applications	Blocking
Target	RPS6KB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214542.
Purity	> 85%
Source	Synthetic
Uniprot ID	P23443, Q8BSK8, P67999
Entrez	6198, 72508, 83840
Expiration Date	12 months from date of receipt.