

Product Datasheet

RPS19 peptide (orb217902)

Description	The peptide is used to block anti-RPS19 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	RPS19
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214546.
Purity	> 85%
Source	Synthetic
Uniprot ID	P39019, Q9CZX8, P17074
Entrez	100910336, 29287, 20085, 6223
Expiration Date	12 months from date of receipt.