

Product Datasheet

CLIP1 peptide (orb217904)

Description	The peptide is used to block anti-CLIP1 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CLIP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214548.
Purity	> 85%
Source	Synthetic
Uniprot ID	P30622, Q922J3
Entrez	6249, 56430
Expiration Date	12 months from date of receipt.