

Product Datasheet

S100A1 peptide (orb217906)

Description	S100A1 peptide
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	S100A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214550.
Purity	> 85%
Source	Synthetic
Uniprot ID	P56565, P23297, P35467
Entrez	20193, 295214, 6271
Expiration Date	12 months from date of receipt.