

Product Datasheet

SCN2B peptide (orb217912)

Description	The peptide is used to block anti-SCN2B antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SCN2B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214556.
Purity	> 85%
Source	Synthetic
Uniprot ID	O60939, P54900, Q56A07
Entrez	72821, 6327, 25349
Expiration Date	12 months from date of receipt.