

Product Datasheet

SCN8A peptide (orb217914)

Description	The peptide is used to block anti-Nav1.6 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SCN8A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214558.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9UQD0 , Q9WTU3 , O88420
Entrez	20273 , 29710 , 6334
Expiration Date	12 months from date of receipt.