

Product Datasheet

MAP2K4 peptide (orb217919)

Description	The peptide is used to block anti-MKK4 antibody reactivity.
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	MAP2K4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214563.
Purity	> 85%
Source	Synthetic
Uniprot ID	P45985, P47809
Entrez	26398, 6416
Expiration Date	12 months from date of receipt.