

Product Datasheet

SFSWAP peptide (orb217923)

Description	The peptide is used to block anti-SWAP antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	SFSWAP
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214567.
Purity	> 85%
Source	Synthetic
Uniprot ID	D3ZTQ1, Q3USH5, Q12872
Entrez	6433, 231769, 304431
Expiration Date	12 months from date of receipt.