

Product Datasheet

SKP2 peptide (orb217932)

Description	The peptide is used to block anti-SKP2 antibody reactivity.
Reactivity	Human, Porcine
Tested Applications	Blocking
Target	SKP2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214576.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q13309
Entrez	6502
Expiration Date	12 months from date of receipt.