

Product Datasheet

SLC2A1 peptide (orb217934)

Description	The peptide is used to block anti-GLUT1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	SLC2A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214578.
Purity	> 85%
Source	Synthetic
Uniprot ID	P17809, P11166, P11167
Entrez	6513, 20525, 24778
Expiration Date	12 months from date of receipt.