

Product Datasheet

GLUT2 peptide (orb217935)

Description	GLUT2 peptide
Reactivity	Human, Rat
Tested Applications	Blocking
Target	SLC2A2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb214579.
Purity	> 85%
Source	Synthetic
Uniprot ID	P12336, P11168
Entrez	6514
Expiration Date	12 months from date of receipt.